



Synaptic  
MEDICAL



# NaviGo<sup>®</sup>

Bi-directional  
Steerable Sheath

[www.synapticmed.com](http://www.synapticmed.com)

# NaviGo®

## Bi-directional Steerable Sheath



- Precise steering control combined with bi-directional deflection for maximum maneuverability
- Ergonomic handle with neutral position indicator
- Advanced valve design provides effective hemostasis
- Excellent sheath-to-dilator transition designed to reduce trauma during transeptal crossing
- Smooth braided shaft provides excellent torqueability, pushability and kink resistance

Catalog Number	Sheath Inner Diameter (F)	Max Guidewire	Curl Type Bi-directional	Usable Length (cm)	Compatible Needle Length (cm)
S08574LS	8.5	.032	Small	71	98
S08574LM	8.5	.032	Medium	71	98
S08574LL	8.5	.032	Large	71	98

Contents: NaviGo steerable sheath, dilator, guidewire



**Synaptic**  
MEDICAL

760.608.8388 | fax 760.608.8528 | sales@synapticmed.com  
[www.synapticmed.com](http://www.synapticmed.com)