



Synaptic  
MEDICAL




**NaviEase<sup>®</sup>**  
Fixed Curve Sheath

[www.synapticmed.com](http://www.synapticmed.com)

# NaviEase<sup>®</sup>

## Fixed Curve Sheath



- Proprietary shaft provides superior torqueability and pushability
- Excellent curve retention improves catheter placement
- Radiopaque marker at tip of sheath enhances fluoroscopic visibility
- Sideport holes on distal end of sheath to ease catheter aspiration and flushing

Catalog Number	Sheath Inner Diameter (F)	Max Guidewire	Curl Type	Usable Length (cm)	Compatible Needle Length (cm)
S0806238R0	8	.038	R0	62	71
S0808138R0	8	.038	R0	81	89
S0808132L1	8	.032	L1	81	89
S0856332L0	8.5	.032	L0	63	71
S0856332L1	8.5	.032	L1	63	71
S0856332L2	8.5	.032	L2	63	71
S0856332LM	8.5	.032	LM	63	71
S0856332LM1	8.5	.032	LM1	63	71
S0856332LM2	8.5	.032	LM2	63	71
S0856332LM3	8.5	.032	LM3	63	71
S0856238R0	8.5	.038	R0	62	71
S0856238R1	8.5	.038	R1	62	71
S0856238R2	8.5	.038	R2	62	71
S0856238RM1	8.5	.038	RM1	62	71
S0856238RM2	8.5	.038	RM2	62	71
S0858132L1	8.5	.032	L1	81	89

Contents: NaviEase sheath, dilator, guidewire



760.608.8388 | fax 760.608.8528 | sales@synapticmed.com  
[www.synapticmed.com](http://www.synapticmed.com)