




Synaptic  
MEDICAL





## Rithm ID™

Fixed Curve  
Diagnostic Catheter

# Rithm ID

## Fixed Curve Diagnostic Catheter



- Designed to optimize torque response for ultimate control
- Provides consistent, long-term wall contact with less repositioning
- Platinum iridium alloy electrodes for improved signal fidelity
- Different sizes and multi curve types

Catalog Number	French Size	Curve Type	Electrodes	Spacing	Usable Length (cm)
F6004005115ART	6	A	4	5-5-5	115
F6004010115ART	6	A	4	10-10-10	115
F6004252115ART	6	A	4	2-5-2	115
F6004252115DRT	6	D	4	2-5-2	115
F6004252115FRT	6	F	4	2-5-2	115
F6010252065PRT	6	P	10	2-5-2	65
F6010282065PRT	6	P	10	2-8-2	65



760.608.8388 | fax 760.608.8528 | sales@synapticmed.com  
[www.synapticmed.com](http://www.synapticmed.com)